
# Welch Allyn Masted Scales

---

---

## Directions for use

Models 5002, 5122, 5702, 6002 (with serial numbers beginning with the letter "U")

© 2018 Welch Allyn. All rights are reserved. To support the intended use of the product described in this publication, the purchaser of the product is permitted to copy this publication, for internal distribution only, from the media provided by Welch Allyn. No other use, reproduction, or distribution of this publication, or any part of it, is permitted without written permission from Welch Allyn.

Welch Allyn assumes no responsibility for any injury to anyone, or for any illegal or improper use of the product, that may result from failure to use this product in accordance with the instructions, cautions, warnings, or statement of intended use published in this manual.

Software in this product is Copyright 2018 Welch Allyn or its vendors. All rights are reserved. The software is protected by United States of America copyright laws and international treaty provisions applicable worldwide. Under such laws, the licensee is entitled to use the copy of the software incorporated with this instrument as intended in the operation of the product in which it is embedded. The software may not be copied, decompiled, reverse-engineered, disassembled, or otherwise reduced to human-perceivable form. This is not a sale of the software or any copy of the software; all right, title, and ownership of the software remain with Welch Allyn or its vendors.

For information about any Welch Allyn product, contact your local Welch Allyn representative: <http://www.welchallyn.com/en/other/contact-us.html>.

# 727271

This manual applies to REF 901109 MASTED SCALES.

DIR 80022302 Ver. C Revision date: 2018-06

Welch Allyn, Inc.  
4341 State Street Road  
Skaneateles Falls, NY 13153 USA

[www.welchallyn.com](http://www.welchallyn.com)

**WelchAllyn**<sup>®</sup>  
Advancing Frontline Care<sup>™</sup>

# Contents

---

<b>Introduction .....</b>	<b>1</b>
Intended use .....	1
Indications for use .....	1
Intended clinical care environments .....	1
Symbols .....	1
About warnings and cautions .....	3
Controls and indicators .....	5
<b>Assembly and setup .....</b>	<b>9</b>
Models 5002 and 5702 stand-on scales .....	9
Model 5122 stand-on scale .....	10
Model 6002 wheelchair scale .....	10
<b>Scale operation .....</b>	<b>11</b>
Power options .....	11
Audible beeps .....	11
Weighing procedure .....	12
Wheelchair weighing procedure (model 6002) .....	12
Calculating BMI .....	13
Print a weight measurement .....	14
Send measurements to a computer .....	14
Send measurements to a Welch Allyn device .....	15
<b>Scale transportation .....</b>	<b>17</b>
Manual transport (models 5002, 5702, and 6002) .....	17
Manual transport (model 5122) .....	17
<b>Storage, cleaning, and maintenance .....</b>	<b>19</b>
Storage .....	19
Cleaning and disinfecting .....	19
Battery replacement .....	19
Change the printer paper .....	20
Maintenance .....	20
General compliance and standards .....	20
Troubleshooting .....	21
<b>Appendices .....</b>	<b>23</b>
Custom setup .....	23
Specifications .....	25
Product configurations .....	27
Approved accessories .....	29
Warranty .....	30
EMC compliance .....	30

# Introduction

---

## Intended use

The Welch Allyn Masted Scales are intended to be used by clinicians for weighing patients up to 660 pounds (300 kg) for model 5122, up to 880 pounds (400 kg) for models 5002 and 6002, and up to 1000 pounds (454 kg) for model 5702.

## Indications for use

The Welch Allyn Masted Scales are used by clinicians to weigh patients. An optional height gauge is available on some models.

Masted scales can make contact with a patient's hands and feet. Contact duration is intended to be limited to less than 30 seconds.

## Intended clinical care environments

The Welch Allyn Masted Scales are intended to be used in the following clinical care environments:

- Hospitals
- Ambulatory care centers
- Physicians' offices
- Other professional medical facilities

## Symbols

### Documentation symbols

**WARNING** The warning statements in this manual identify conditions or practices that could lead to illness, injury, or death.

**Note**

Warning symbols will appear with a grey background in a black and white document.

**Caution** The caution statements in this manual identify conditions or practices that could result in damage to the equipment or other property, or loss of data.

Mandatory - Consult Directions for Use

## Shipping, storing, and environment symbols

Temperature limitation

Keep dry

Separate collection of Electrical and Electronic Equipment. Do not dispose as unsorted municipal waste.

Fragile

Humidity limitation

Atmospheric pressure limitation

This way up

Stacking limit by number

Recyclable

Atmospheric pressure limitation

## Miscellaneous symbols

Manufacturer

Reorder Number

Product Identifier

For indoor use only

Serial Number

**R<sub>x</sub> ONLY**

Prescription only or "For Use by or on the order of a licensed medical professional"

Battery

Direct current (DC)

USB

On/Off Pushbutton

	Class II equipment		Type B applied part <b>Note</b> The entire scale is considered an applied part.
	Rated power input, DC		Mass
	With respect to electrical shock, fire, and mechanical hazards only In accordance with: ANSI/AAMI ES60601-1:2005/(R)2012 and A1:2012, C1:2009/(R)2012 AND A2:2010/(R)2012, CAN/CSA C22.2 No. 60601-1:14 IEC 60601-1 Ed. 3.1		Global Trade Item Number

## About warnings and cautions

Caution statements can appear on the Welch Allyn Masted Scales, the packaging, the shipping container, or in this *Directions for use*.

The Welch Allyn Masted Scales are safe for patients and clinicians when used in accordance with the instructions and caution statements presented in this *Directions for use*.

Before using the device, you must familiarize yourself with all cautions, with the steps to power up the device, and with the sections of this *Directions for use* that pertain to your use of the device. In addition to reviewing the general cautions presented in the next section, you must also review the more specific cautions that appear throughout the manual in conjunction with setup/startup, operation, and maintenance tasks. No additional training is required.

- Failure to understand and observe any warning statement in this manual could lead to patient injury or illness.
- Failure to understand and observe any caution statement in this manual could lead to damage to the equipment or other property, or loss of patient data.

## General warnings and cautions

**WARNING** Patient or operator injury risk. Maintain control of the patient at all times, including when the patient is stepping on and off the scale. Failure to do so could result in a fall risk or other serious injury to you or the patient.

**WARNING** Patient or operator injury risk. Make sure the scale is positioned on a stable surface and in a way to avoid shock or bumping hazards.

**WARNING** Patient injury risk. Make sure the scale is set in the proper measurement units for your facility. Make sure you record the measurement units as they are displayed.

**WARNING** Electric shock hazard. Use only a Welch-Allyn approved power supply. The use of an unapproved power supply could increase chassis or patient leakage currents.

**WARNING** Electric shock hazard. Use only a Welch-Allyn approved power supply. The use of an unapproved power supply could cause electric harm and shock to you or the patient.

**WARNING** Electric shock hazard. All signal input and output (I/O) connectors are intended for connection of only devices complying with IEC 60601-1, or other IEC standards (for example, IEC 60950), as applicable to the scale. Connecting additional devices to the scale could increase chassis or patient leakage currents. To maintain operator and patient safety, consider the requirements of IEC 60601-1-1. Measure the leakage currents to confirm that no electric shock hazard exists.

**WARNING** Electric shock hazard. Make sure that the USB (SIP/SOP) port and the patient are never touched or come in contact at the same time.

**WARNING** Patient or operator injury risk. No modification of this equipment is allowed.

**CAUTION** Make sure that you routinely perform general maintenance and equipment safety checks on your scale. Remove the scale from service when you notice damaged power cable.

**CAUTION** Do not use this scale to transport patients or items.


**CAUTION** The scale may power down due to electrostatic discharge. There is no risk associated with this action. Turn the scale on and it will function normally.


# Controls and indicators

## Scale


### Models 5002 and 5702

- 
- 1 Steering handle/hand grip
  - 2 Control panel
  - 3 Mast
  - 4 Hand rail
  - 5 Wheels
  - 6 Weighing platform
- N/A USB port (located on the side of the mast)
- N/A Printer (optional, located on the front of the mast)
- N/A Battery access door (located on the rear of the mast)
- N/A Power input (located on the sides of the mast)
- N/A Height gauge (optional)
- N/A Wall-mounted height gauge (in and cm) (optional)
-

**Model 6002**

- 
- 1 Steering handle/hand grip
  - 2 Control panel
  - 3 Mast
  - 4 Hand rail (optional)
  - 5 Wheels
  - 6 Weighing platform
  - 7 Wheelchair ramp
- N/A USB port (located on the side of the mast)
  - N/A Printer (optional, located on the front of the mast)
  - N/A Battery access door (located on the rear of the mast)
  - N/A Power input (located on the sides of the mast)
  - N/A Height gauge (optional)
  - N/A Wall-mounted height gauge (in and cm) (optional)
-

**Model 5122**

---

1 Control panel

2 Mast (battery access door located on front)

3 Weighing platform

N/A USB port (optional, located on the side of the mast)


N/A Power input (located on the sides of the mast)

N/A Wall-mounted height gauge (in and cm) (optional)

N/A Carrying handle (located on the rear of the mast)

---

## Control panel buttons

Button	Function
ZERO	Press this button to reset the scale value to zero. This button can also be used for a tare function to tare out the weight of a chair, wheelchair, or other components you wish to be deducted from the patient's weight
Green LED indicator	This area displays the patient measurement, as well as indicates the measurement units, AC power, low battery, prior weight, and data transmission.
Power button 	Press this button to turn the scale on. Press and hold this button to turn the scale off.
SEND DATA	Press this button to send measurement data when the scale is connected to a computer, a Welch Allyn Connex Vital Signs® Monitor or a Welch Allyn Connex® Spot Monitor.
PRINT	Press this button to print the current measurements.
LB / KG	Press this button to change the scale between kilograms and pounds.
REWEIGH	Press this button to retake the patient's weight or to perform a new weighing cycle.
HEIGHT / BMI	Press this button to calculate a BMI measurement or enter a patient's height.
RECALL	Press this button to recall the previously stored scale reading.

# Assembly and setup

---

The scale is fully assembled at the factory and tested. In some instances, the scales are shipped assembled. However, due to size and shipping restrictions, most scales are partially disassembled. If assembly is required, it can be done with common hand tools. When special tools are required, they are supplied with the scale. Ensure that all bolts and fasteners are securely tightened.

Use the following instructions to assemble and setup your scale for use.

**Note** If your scale requires assembly, make sure that the serial number of the mast and the serial number of the base match when assembling.

## Models 5002 and 5702 stand-on scales

1. Carefully remove all parts from the shipping carton and inspect the components for shipping damage.
2. Place the scale base on the floor or on a workbench.
3. Remove the bolts and washers from the rear of the scale base.
4. Insert the mast into the base. Ensure that the readout display is facing the patient when the patient stands on the scale.
5. Fasten the mast with the bolts and washers removed in Step 3.
6. Tighten the bolts securely to ensure that a solid connection is made.

## Weighing deck mounting (model 5702)

1. Position the weighing deck so that the two sets of holes are closest to the mast.
2. Line the weighing deck up with the mounting holes in the aluminum plate.
3. Install the four flat head screws through the weighing deck.
4. Securely tighten the supplied mounting screws with the provided hex key.

## Handrail mounting

1. Attach the handrail to the top of the weighing platform using the button head screws and the provided hex key.  
Three screws are required for each mounting plate (six total screws).
2. Tighten all six screws securely to ensure a safe and stable handrail.

## Height gauge assembly (purchasable option)

If your scale was ordered with a height gauge, follow the installation instructions included with the height gauge.

## Model 5122 stand-on scale

The 5122 model scale is shipped assembled, so no assembly is required. It operates the same as the 5002 model. Handrails are not available for this scale.

## Model 6002 wheelchair scale

1. Carefully remove all parts from the shipping carton and inspect the components for shipping damage.
2. Place the scale base on the floor or on a workbench.
3. Remove the bolts and washers from the rear of the scale base.
4. Insert the mast into the base. Ensure that the readout display is facing the patient when the patient stands on the scale.
5. Fasten the mast with the bolts and washers removed in Step 3.
6. Tighten the bolts securely to ensure a solid connection is made.

## Weighing deck mounting

1. Position the weighing deck so that the two sets of holes under the flaps are closest to the mast.
2. Line the weighing deck up with the mounting holes in the aluminum plate.
3. Install the four flat head screws through the weighing deck.
4. Securely tighten the supplied mounting screws with the provided hex key.

## Handrail mounting (purchasable option)

If your scale was purchased with an optional handrail, follow the installation instructions included with the handrail.

## Height gauge assembly (purchasable option)

If your scale was ordered with a height gauge, follow the installation instructions included with the height gauge.

# Scale operation

---

## Power options

The scale can be used with battery power or with AC power.

Press the power button to turn on the scale. When connected to AC power, the green AC power light illuminates.

### Battery power

**CAUTION** Use only size D disposable alkaline batteries. The use of any other battery will void the warranty.

Press  to turn on the scale.

The scale automatically turns off after 45 seconds of non-use when operating on battery power.

The **LOW BATTERY** indicator illuminates when the battery power is low.

**bAttrY** displays when the batteries are unable to properly power the scale.

Replace the batteries when either of these messages are displayed.

### AC power

When plugged in, the scale continuously runs on AC power automatically, depending on your settings.

To disconnect the scale from AC power, either remove the power adaptor from the wall or the side of the scale. Make sure to position the scale so that it can be easily disconnected from AC power.

The scale switches to battery power when AC power is not present.

## Audible beeps

The following beeps occur when using the scale.

**Note** These beeps do not occur if the beeper option has been disabled in the setup menu.

---

Sound	Meaning
-------	---------

One beep	Acknowledgment of a button press.
Two beeps	The scale is zeroed and ready to take a measurement.
Three beeps	The weight reading is obtained and displayed.
No print, four short beeps	A print was requested, but valid weight is not displayed on the scale, or <b>PRINT</b> was pressed while the printer is in use.
No print, four long beeps	There is a problem with the printer. Contact your service department.
Four long beeps	The battery is too weak to operate the scale. The scale has entered setup mode.
Continuous long beep	There is a problem with the scale. See the troubleshooting section.

---

## Weighing procedure

Make sure to select a consistent time of day to obtain a weight measurement.

Position the scale in a convenient and stable location.

1. Press **ON** to turn on the scale.  
The scale sounds one beep.
2. Wait for the scale to display **0.0**.  
The scale sounds two beeps.
3. Place the patient on the scale and make sure that his or her weight is evenly distributed.  
The scale sounds three beeps and the patient's weight is displayed.
4. If necessary, press **REWEIGH** to perform another weight measurement.
5. If necessary, press **RECALL** to display the previously stored scale measurement.  
The **PRIOR WEIGHT** LED flashes.

**Note** A measurement is stored in the scale until a new weight is taken and zeroes out.

A series of moving dashes indicates that a tared weight is stored in the scale memory, or that there is a negative weight on the scale. Press **ZERO** to clear the tare value and reset the scale.

If **O-LOAD** is displayed, the capacity of the scale has been exceeded.

## Wheelchair weighing procedure (model 6002)

1. Make sure that the scale ramps are down.
2. Press the power button to turn on the scale.  
The scale sounds one beep.
3. Wait for the scale to display **0.0**.  
The scale sounds two beeps.

4. Place the empty wheelchair on the scale and make sure the weight is evenly distributed.  
A series of dashes moves across the display. The scale sounds three beeps and the empty wheelchair weight is displayed.
5. Press **ZERO** to tare out the weight of the wheelchair.  
A series of dashes moves across the display.
6. When the scale displays **0.0**, remove the wheelchair.  
A series of dashes moves across the display. This means a tared weight is stored in the scale memory.
7. Place the patient in the wheelchair and place both on the scale. Make sure the weight is evenly distributed.  
The scale sounds three beeps and the patient's weight is displayed.
8. If necessary, press **REWEIGH** to perform another weight measurement with the current tare value.
9. To remove the tare value and reset the scale, remove all weight from the scale and press **ZERO**.
10. If necessary, press **RECALL** to display the previously stored scale measurement.  
The **PRIOR WEIGHT** LED flashes.

**Note** A measurement is stored in the scale until a new weight is taken and zeroes out.

A series of moving dashes indicates that a tared weight is stored in the scale memory, or that there is a negative weight on the scale. Press **ZERO** to clear the tare value and reset the scale.

If **O-LOAD** displays, the capacity of the scale has been exceeded.

If **OFF SCALE PLEASE** – displays, remove the patient from the scale. Wait for the scale to display **0.0** and retry the measurement.

**Note** When using battery power, the scale does not save tare values after powering down. See the Custom setup section to change the number of seconds before the scale turns off when operating on battery power.

## Calculating BMI

Use the following procedure to calculate a patient's body mass index (BMI).

1. Press the power button to power on the scale.  
The scale sounds one beep.
2. Wait for the scale to display **0.0**.  
The scale sounds two beeps.
3. Have the patient step on the scale, and make sure his or her weight is evenly distributed.  
The scale sounds three beeps and the patient's weight is displayed.
4. If necessary, press **LB/KG** to change the unit of measure on the display.

**Note** If you purchased the kilogram-only option, the scale displays Lb OFF when **LB/KG** is pressed.

5. Press **HEIGHT/BMI**.

A default value is displayed, and the height units LED is illuminated.

6. Use the height gauge to measure the patient's height.
7. Use the up and down arrows to enter the patient's height. Press and hold an arrow up or down to raise or lower the patient's height value at an accelerated rate.

**Note** If the scale is set in kilograms-only mode, the height is displayed in centimeters. If the scale is set in pounds-only or in kilograms and pounds mode, the height is displayed in inches.

8. Press **HEIGHT/BMI**.

The scale displays the BMI measurement and the BMI unit LED is illuminated. Press **HEIGHT/BMI** to cycle through the weight, height, and BMI values.

**Note** The scale does not calculate BMI if the patient's height or weight are not captured.

## Print a weight measurement

This section applies to scales configured with an optional printer at time of purchase. If the printer is configured to Auto Print (see PrtOpt\* in the Custom Setup section), a printed record is automatically generated after a measurement is taken.

1. Take a weight measurement as described.  
A single beep sounds.
2. Press **PRINT**.  
A single beep sounds and the measurement is printed.
3. To print a previously stored scale measurement, press **RECALL** and then press **PRINT**.

The following errors can occur:

- Four short beeps indicate an invalid print request. Make sure that the weight is properly displayed on the scale, or that a measurement is not already in the printing process.
- Four long beeps indicate an internal problem with the printer. See the Troubleshooting section.

## Send measurements to a computer

The scale can connect to a computer using the Microsoft® Windows® 10 operating system.

Detailed instructions for this feature are available from Welch Allyn. Go to <http://www.welchallyn.com/en/other/contact-us.html> to find your local representative.

**CAUTION** Any computer connected to the scale must be running on battery power, a 60601-1-compliant power supply, or a 60601-1-compliant isolation transformer.

1. Connect the USB Type B male side of the cable to the USB port on the side of the mast.
2. Connect the USB Type A male side of the cable into a USB port on your computer.  
Windows automatically installs the necessary driver.
3. Start the customer-supplied data collection software configured to receive data on the COM port assigned by Windows.
4. Press **SEND DATA** on the scale to send the current measurement to the computer.

## Send measurements to a Welch Allyn device

The scale can connect to a Welch Allyn Connex<sup>®</sup> Vital Signs Monitor or a Welch Allyn Connex<sup>®</sup> Spot Monitor.

**CAUTION** Any device connected to the scale must be running on battery power, a 60601-1-compliant power supply, or a 60601-1-compliant isolation transformer.

1. Connect the USB Type B male side of the cable to the USB port on the side of the mast.
2. Connect the USB Type A male side of the cable into a USB port on the rear of your device.
3. Press **SEND DATA** on the scale to send the current measurement to the device.

# Scale transportation

---

## Manual transport (models 5002, 5702, and 6002)

1. Unplug the power adapter cable from the scale.
2. Tip the scale backward to engage the wheels on the rear of the scale.
3. Use the only steering handle to push the scale to the new location.
4. At the new scale location, do the following:
  - a. Make sure the scale is positioned on a flat surface and that all four corners of the weighing platform touch the ground. This acts as the scale's braking system.
  - b. Make sure to plug the power adapter cable into the wall and into the side of the scale.
  - c. Make sure the power adapter cable does not pose a tripping hazard and is easily accessible.

## Manual transport (model 5122)

1. Unplug the power adapter cable from the scale.
2. Use the carrying handle located on the rear of the mast to lift the scale and carry it to the new location.
3. At the new scale location, do the following:
  - a. Make sure the scale is positioned on a flat surface and that all four corners of the weighing platform touch the ground. This acts as the scale's braking system.
  - b. Make sure to plug the power adapter cable into the wall and into the side of the scale.
  - c. Make sure the power adapter cable does not pose a tripping hazard and is easily accessible.

# Storage, cleaning, and maintenance

---

## Storage

The scale should be stored in a convenient storage facility or closet. Make sure to not leave the scale in an area where someone could trip on it or run into it.

## Cleaning and disinfecting

**Note** Any visible soil must be removed before cleaning or disinfecting.

1. Use one of the following products when cleaning your scale:
  - PDI Super Sani-Cloth® (disinfectant) and Sani-Cloth® Plus (disinfectant)
  - Metrex CaviWipes™ (disinfectant)
  - A clean cloth moistened with 70 percent isopropyl alcohol
  - Clorox Healthcare® Bleach Germicidal Wipes (disinfectant) or Sani-Cloth® Bleach Wipes (disinfectant)
2. For optimal results, follow the manufacturer's instructions and your facility's procedure to clean the scale. Make sure to properly dispose of any used cleaning supplies.

## Battery replacement


**WARNING** Electronic shock hazard. Do not replace the battery under the following conditions:

- in the vicinity of the patient
- with wet hands
- in the presence of flammable anesthetics

**CAUTION** Use only size D disposable alkaline batteries. The use of any other battery will void the warranty.

1. Make sure the scale is powered down.
2. Using a Phillips screwdriver, remove the two fasteners on the access door and remove the door.
3. Remove the old batteries and install six new size D batteries in the battery holder. Make sure to follow the polarity instructions.

4. Re-attach the access door by installing the two screws.

## Change the printer paper

The printer is a purchasable option. Perform the following procedure if your scale is equipped with a printer.

1. Use the finger hold to lift the clear window on the front of the printer to open the printer door.
2. Remove the old printer paper spindle.
3. Replace the printer paper. Make sure the new printer paper is advancing from the top of the roll.
4. Pull a small piece of the printer paper through the top of the printer door, and close the printer door.

## Maintenance

Perform the following preventive maintenance quarterly to keep your scale in working order, or as required based on usage or per your facility maintenance schedule, whichever comes first.

1. Have your service department check the calibration annually or as required. Calibration instructions are provided in the scale service manual.
2. Inspect the scale for cracks or loose mounting hardware. Replace or repair as necessary.
3. Visually inspect the scale enclosure for damage or loose or missing hardware. Replace or repair as necessary.
4. Inspect the power adapter cable for abrasions or other signs of wear.
5. Do not expose the scale to excessive water or moisture.
6. Do not store the scale where heavy objects can be placed on it.
7. Replace the batteries annually or as required.
8. When storing the scale, remove the batteries from the battery area. Batteries can corrode over a period of time. Make sure to check the batteries before putting the scale back into use.
9. Do not service or perform maintenance while the scale is in use with a patient.

## General compliance and standards

**Directive 2002/96/EC-WEEE:**  
**Disposal of noncontaminated electrical and electronic equipment**

This product and its accessories must be disposed of according to local laws and regulations. Do not dispose of this product as unsorted municipal waste. Prepare this product for reuse or separate collection as specified by Directive 2002/96/EC of the European Parliament and the Council of the European Union on Waste Electronic and Electrical Equipment (WEEE). If this product is contaminated, this directive does not apply.

For more specific disposal or compliance information, see [www.welchallyn.com/weee](http://www.welchallyn.com/weee), or contact Welch Allyn Customer Service at +44 207 365 6780.

## Troubleshooting

This section presents tables of information to help you troubleshoot issues with the scale.

To use these tables, locate the specific problem with your scale in the left column of the table. The remainder of the row explains possible causes and suggests actions that can resolve the issue.

Issue	Possible cause	Suggested action
The weight reading is not accurate	There may be a mechanical obstruction on the platform	Make sure the platform is not touching a foreign object or something restricting its downward motion.
	The cable between the platform and the instrument circuit board may not be connected	Connect the cable.
	The cable on the mast may not be connected	Connect the cable.
The weight reading takes too long to display	The platform may be in motion	Wait for the platform to settle.
The <b>LOW BATTERY</b> indicator is illuminated	The scale batteries are low on power	Replace the batteries.
The scale displays <b>LoBATT</b>	The batteries do not have enough power to print	Replace the batteries.
The scale displays <b>bAttrY</b>	The batteries do not have enough power to operate the scale	Replace the batteries.
The scale displays <b>CABLE</b>	There may be a cable failure	Contact your service department.
The scale displays <b>O-LOAD</b>	The object on the scale may be larger than the maximum value allowed	If the weight value is within the specified range, contact your service department.
The reading does not change when weight is applied to the scale	The load cell transducer, connector, or cable may be defective	Contact your service department.
The scale buttons are not responding	The buttons may be damaged	Check for visible signs of damage.
		Contact your service department.

Issue	Possible cause	Suggested action
The error message <b>ERR FD</b> is displayed	The scale's factory database has failed	Contact your service department.
The error message <b>ERR DB</b> is displayed	The scale's application database has failed	Contact your service department.
The error message <b>ERRoR</b> is displayed	The scale's flash memory has failed	Contact your service department.
The error message <b>ERR AD</b> is displayed	The measurement board cannot be read	Contact your service department.
The error message <b>ERR AV</b> is displayed	The measurement board is the wrong version	Contact your service department.
The error message <b>ERR CL</b> is displayed	Measurement board error, the weight display is blocked	Contact your service department.
The printed paper tears poorly or unevenly	The tear bar may be damaged	Check for signs of damage such as missing teeth or cracks. If damaged, contact your service department.

# Appendices

---

## Custom setup

You can customize the scale to best suit your needs. Configurations that can be selected include the automatic shut-off time, weighing units, weight resolution, and more.

### Enter custom setup mode

1. Make sure that the scale is powered down.
2. Press and hold **ZERO** and **REWEIGH** while pressing the power button.  
**0000** is displayed.
3. Enter the PIN **9821** using the following keys:
  - Press **REWEIGH** to adjust the selected value up or press **RECALL** to adjust the selected value down.
  - Press **LB/KG** to confirm the entry and move to the next value.
4. Once the correct PIN is entered, press **LB/KG** to enter the setup mode  
The scale displays **Set-UP** and produce a series of four long beeps.

### Set the options

Enter the custom setup as described.

1. Press **LB/KG** to scroll through the custom setup options.  
The options displayed with each press are as follows:  
Press **REWEIGH** to adjust these options up or press **RECALL** to adjust these options down. Press **LB/KG** to accept the current value and move on.
2. Press and hold the power button to power down the scale. Power the scale back on to take measurements with the new values.

**Note** Options indicated with an asterisk (\*) require an additional press of **LB/KG** to change the value.

Option displayed	Feature
SOft	This option displays the software version of your scale.
dAtE	This option displays the release date of the software. The format is MM.DD.YY.

SCALE	This option displays the model number of the scale.
AutOFF*	This displays the number of seconds before the scale turns off when operating on battery power. The <b>Cont</b> value prevents the scale from turning off automatically.
ACCONt*	This option causes the scale to remain on when plugged in to AC power.
rES	Pressing <b>REWEIGH</b> or <b>RECALL</b> switches between the following options: 0.1 lb / 0.05 kg 0.1 lb / 0.1 kg 0.02 lb / 0.1 kg 0.05 lb / 0.2 kg 1.0 lb / 0.5 kg  <b>Note</b> If you change the resolution, restart the scale and press <b>ZERO</b> to reset the scale value to zero.
UnitS*	This option allows you to change the displayed weight unit.  <b>Note</b> Do not change the scale units if you have purchased the kilogram-only option. KILoS* This option allows the weight to display in kilograms. POUNDS* This option allows the weight to display in pounds.
bEEPEr*	This option determines the audible signal that occurs when a front panel button is pressed.
rEcALL	This option allows you to turn on the scale in Recall mode. If this option is set to OFF, the <b>RECALL</b> button does not operate.
USb*	This option allows you to send weight values to a computer or a Welch Allyn device.
prNtEr*	This option allows you to turn the printer on or off.  <b>Note</b> The printer is a purchasable option. This option is only applicable to scales equipped with a printer.
Prtunt*	This option determines whether pounds or kilograms are displayed on a printout. Press <b>REWEIGH</b> or <b>RECALL</b> to switch between the following options: If the pounds and kilograms indicators are illuminated on the control panel, the printed output is determined by whichever unit is selected on the control panel. If the kilogram indicator is illuminated, it will only print in kilograms, regardless of the unit selected on the control panel. If the pounds indicator is illuminated, it will only print in pounds, regardless of the unit selected on the control panel.
PrtOPt*	Press <b>REWEIGH</b> or <b>RECALL</b> to switch between the following options: <b>PnIPrt</b> : The current weight is printed when PRINT is pressed on the control panel. <b>AutPrt</b> : The current weight is automatically printed when a weight reading occurs.
SEt-UP	This option is displayed when you have cycled through all the options. Press and hold the power button to power down the scale.

# Specifications

This device is powered by an external power supply. When used together, the device and the power supply are considered ME equipment.

Item	Specification
Scales	5002, 5122, 5702, 6002
Accuracy	0.4 lb for weights between 4 lb to 40 lb 1% for weights between 40 lb to the weighing capacity of the respective scale models
Resolution (factory default)	0.1 lb / 100 gm
Environmental specifications	<p><b>Operation</b></p> <ul style="list-style-type: none"> <li>• Temperature : 10 °C to 40 °C (50 °F to 104 °F)</li> <li>• Humidity: 10% RH to 90% RH, non-condensing</li> <li>• Atmospheric pressure: 700 hPa to 1060 hPa</li> </ul> <p><b>Shipping</b></p> <ul style="list-style-type: none"> <li>• Temperature: -20 °C to 50 °C (-4 °F to 122 °F)</li> <li>• Humidity: 15% to 95% RH</li> <li>• Atmospheric pressure: 700 hPa to 1060 hPa</li> </ul> <p><b>Storage</b></p> <ul style="list-style-type: none"> <li>• Temperature -35 °C, humidity uncontrol</li> <li>• Temperature 30 °C, humidity 90%</li> <li>• Temperature 55 °C, humidity 30%</li> </ul>
Protection rating	IPX0
Weighing range	Models 5002 and 6002: 4 lb to 880 lb (1.8 kg to 400 kg) Model 5702: 4 lb to 1000 lb (1.8 kg to 454 kg) Model 5122: 4 lb to 660 lb (1.8 kg to 300 kg)
Readout	Large, bright light-emitting diode display. Selectable weight displayed in pounds or kilograms. Indicator light to show selected units.
Power source	<p>Battery</p> <ul style="list-style-type: none"> <li>• 6 D alkaline batteries (1.5V x 6)</li> </ul> <p>Power adapter:</p> <ul style="list-style-type: none"> <li>• Medical grade UL listed external power supply</li> <li>• Nominal input voltage: 100-240V AC <math>\pm 10\%</math></li> <li>• Nominal input frequency: 50-60 Hz</li> <li>• Nominal input current: 0.4-0.2A<sub>rms</sub> @ max load</li> <li>• Standby power consumption at <math>U_{in}</math>: 230V AC : <math>\leq 0.1W</math></li> <li>• Nominal output volt age: <math>U_{out}</math> : 12V DC +5% / -5% <math>U_{Br}</math> : <math>\leq 150mVpp</math></li> <li>• Nominal output current: <math>I_{out}</math> : 1500mA</li> </ul>
Platform dimensions	5002: 18.25 in. x 19.7 in. x 2.3 in.

Item	Specification
	5702: 24 in. x 26 in. x 2.5 in. 6002: 24 in. x 28 in. x 2.5 in. 5122: 16 in. x 14 in. x 1.6 in.
Scale weight	5002: 62 lb (28.1 kg) 5122: 22 lb (9.9 kg) 5702: 76 lb (34.5 kg) 6002: 78 lb (35.4 kg)

## Technical description

The Welch Allyn Masted Scales use load cell transducers to convert analog weight applied to the weighing platform into a digital, electrical signal for display on the readout.

Welch Allyn will provide, on request, circuit diagrams, component part lists, descriptions, calibration instructions, and other information to authorized service personnel.

## Product configurations

The Welch Allyn Masted Scales are available in the following configurations:

### Model 5002, Mobile stand-on scale

Model	Description
5002-KP-B	Mobile stand-on scale, kg only, printer, dual power (battery and AC)
5002-KX-B	Mobile stand-on scale, kg only, dual power (battery and AC)
5002-KX-X	Mobile stand-on scale, kg only, battery powered
5002-XP-B	Mobile stand-on scale, user-selectable units (lb/kg), printer, dual power (battery and AC)
5002-XX-B	Mobile stand-on scale, user-selectable Units (lb/kg), dual power (battery and AC)
5002-XX-X	Mobile stand-on scale, user-selectable units (lb/kg), battery powered

### Model 5122, Low-profile stand-on scale

Model	Description
5122-K-X	Low-profile stand-on scale, kg only, battery powered
5122-X-X	Low-profile stand-on scale, user-selectable units (lb/kg), battery powered

### Model 5702, Mobile bariatric stand-on scale

Model	Description
5702-KP-B	Mobile bariatric stand-on scale, kg only, printer, dual power (battery and AC)
5702-KX-B	Mobile bariatric stand-on scale, kg only, dual power (battery and AC)
5702-KX-X	Mobile bariatric stand-on scale, kg only, battery powered
5702-XP-B	Mobile bariatric stand-on scale, user-selectable units (lb/kg), printer, dual power (battery and AC)
5702-XX-B	Mobile bariatric stand-on scale, user-selectable units (lb/kg), dual power (battery and AC)
5702-XX-X	Mobile bariatric stand-on scale, user-selectable units (lb/kg), battery powered

### Model 6002, Wheelchair scale

<b>Model</b>	<b>Description</b>
6002-KP-B	Wheelchair scale, kg only, printer, dual power (battery and AC)
6002-KX-B	Wheelchair scale, kg only, dual power (battery and AC)
6002-KX-X	Wheelchair scale, kg only, battery powered
6002-XP-B	Wheelchair scale, user-selectable units (lb/kg), printer, dual power (battery and AC)
6002-XX-B	Wheelchair scale, user-selectable units (lb/kg), dual power (battery and AC)
6002-XX-X	Wheelchair scale, user-selectable units (lb/kg), battery powered

## Approved accessories

### Model 5002

Part number	Description
845010	Height gauge (in and cm)
845010M	Height gauge (metric only)
845010W	Wall-mounted height gauge (in and cm)
412713	Thermal printer paper (box, 15 rolls)
726551	Power supply

### Model 5122

Part number	Description
845010W	Wall-mounted height gauge (in and cm)
726551	Power supply

### Model 5702

Part number	Description
845010B	Height gauge for 5702 (in and cm)
845010W	Wall-mounted height gauge (in and cm)
412713	Thermal printer paper (box, 15 rolls)
726551	Power supply

### Model 6002

Part number	Description
60224	Handrail
845010	Height gauge (in and cm)
845010M	Height gauge (metric only)
845010W	Wall-mounted height gauge (in and cm)
412713	Thermal printer paper (box, 15 rolls)
726551	Power supply

## Warranty

Welch Allyn will warranty the weight scale to be free of defects in material and workmanship and to perform in accordance with manufacturer specifications for the period of one year from the date of retail purchase.

The warranty period shall start on the date of purchase. The date of purchase is: 1) the invoiced ship date if the device was purchased directly from Welch Allyn, 2) the date specified during product registration, 3) the date of purchase of the product from a Welch Allyn authorized distributor as documented from a receipt from said distributor.

This warranty does NOT cover damages caused by misuse or abuse, including but not limited to:

- Failure caused by unauthorized repairs or modifications
- Damage caused by shock or dropping during transportation
- Damage caused by improper use of the power supply
- Failure caused by improper operation not consistent with the instructions stated in this *Directions for use*

Should this device require maintenance (or replacement at our option) under warranty, contact your local Welch Allyn representative: <http://www.welchallyn.com/en/other/contact-us.html>

## EMC compliance

Special precautions concerning electromagnetic compatibility (EMC) must be taken for all medical electrical equipment. This device complies with IEC 60601-1-2 Ed 4.0 / EN 60601-1-2:2015..

- All medical electrical equipment must be installed and put into service in accordance with the EMC information provided in this document, the *Welch Allyn® Masted Scales Directions for use*.
- Portable and mobile RF communications equipment can affect the behavior of medical electrical equipment.

The device complies with all applicable and required standards for electromagnetic interference.

- It does not normally affect nearby equipment and devices.
- It is not normally affected by nearby equipment and devices.
- It is not safe to operate the scale in the presence of high-frequency surgical equipment.
- However, it is good practice to avoid using the device in extremely close proximity to other equipment.

## Emissions and immunity information

---

### Electromagnetic emissions

---

The Welch Allyn Masted Scale is intended for use in the electromagnetic environment specified below. The customer or user of the Welch Allyn Masted Scale should assure that it is used in such an environment.

---

**Electromagnetic emissions**

<b>Emissions test</b>	<b>Compliance</b>	<b>Electromagnetic environment - guidance</b>
Radiated emissions CISPR 11	Group 1 /Class B	The Welch Allyn Masted Scale uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
Conducted emissions CISPR 11	Class B	The Welch Allyn Masted Scale is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

**Electromagnetic immunity**

The Welch Allyn Masted Scale is intended for use in the electromagnetic environment specified below. The customer or the user of the Welch Allyn Masted Scale should assure that it is used in such an environment.

<b>Immunity test</b>	<b>IEC 60601 test level</b>	<b>Compliance level</b>	<b>Electromagnetic environment - guidance</b>
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode	±1 kV differential mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<0 % $U_T$ (>100 % dip in $U_T$ ) for 1 cycle @ 0° 70 % $U_T$ (30 % dip in $U_T$ ) for 25/30 cycles @ 0° <0 % $U_T$ (>100 % drop in $U_T$ ) for 0.5 cycles @ 0°, 45°, 90°, 135°, 180°, 225°, 270°, & 315° 0% for 250/300 cycles @ 0°	<0 % $U_T$ (>100 % dip in $U_T$ ) for 1 cycle @ 0° 70 % $U_T$ (30 % dip in $U_T$ ) for 25/30 cycles @ 0° <0 % $U_T$ (>100 % drop in $U_T$ ) for 0.5 cycles @ 0°, 45°, 90°, 135°, 180°, 225°, 270°, & 315° 0% for 250/300 cycles @ 0°	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Welch Allyn Masted Scale requires continued operation during power mains interruptions, it is recommended that the Welch Allyn Masted Scale be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

**Electromagnetic immunity**

Note:  $U_T$  is the a.c. mains voltage prior to application of the test level.

**Electromagnetic immunity**

The Welch Allyn Masted Scale is intended for use in the electromagnetic environment specified below. The customer or the user of the Welch Allyn Masted Scale should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
			Portable and mobile RF communications equipment should be used no closer to any part of the Welch Allyn Masted Scale, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
			<b>Recommended separation distance</b>
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 V	$d = (1.17) \sqrt{P}$
	6Vrms ISM Band frequency (MHz) : 6.675 – 6.795, 13.553 – 13.567, 26.957 – 27.284 & 40.66 – 40.70	6V	
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz	3 V/m	$d = (1.17) \sqrt{P}$ 80 to 800 MHz
			$d = (2.33) \sqrt{P}$ 800 MHz to 2.7 GHz

where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and  $d$  is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey<sup>a</sup>, should be less than the compliance level in each frequency range<sup>b</sup>. Interference may occur in the vicinity of equipment marked with the following symbol:

**Table-9 of IEC60601-1-2, 4th Ed.**

27 V/m 380 MHz to 390 MHz	27 V/m 380 MHz to 390 MHz
28 V/m 430 MHz to 470 MHz	28 V/m 430 MHz to 470 MHz

---

**Electromagnetic immunity**

---

800 MHz to 960 MHz	800 MHz to 960 MHz
1700 MHz to 1990 MHz	1700 MHz to 1990 MHz
2400 MHz to 2570 MHz	2400 MHz to 2570 MHz
9 V/m	9 V/m
704 MHz to 787 MHz	704 MHz to 787 MHz
5100 MHz to 5800 MHz	5100 MHz to 5800 MHz

---

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Welch Allyn Masted Scale is used exceeds the applicable RF compliance level above, the Welch Allyn Masted Scale should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Welch Allyn Masted Scale.

<sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

---

---

**Recommended separation distances between portable and mobile RF communications equipment and the Welch Allyn Masted Scale**

---

The Welch Allyn Masted Scale is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the Welch Allyn Masted Scale can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Welch Allyn Masted Scale as recommended below, according to the maximum output power of the communications equipment.

---

**Separation distance according to frequency of transmitter (m)**

---

Rated max. output power of transmitter (W)	150 kHz to 80 MHz $d = (1.17) \sqrt{P}$	80 MHz to 800 MHz $d = (1.17) \sqrt{P}$	800 MHz to 2.7 GHz $d = (2.33) \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.20	1.20	2.30
10	3.70	3.70	7.40
100	12	12	23

---

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

---

---

**Recommended separation distances between portable and mobile RF communications equipment and the Welch Allyn Masted Scale**


---

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

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