
Welch Allyn Pediatric/Infant Scale

Directions for use

Model 4802D (with serial numbers beginning with the letter "U")

© 2019 Welch Allyn. All rights are reserved. To support the intended use of the product described in this publication, the purchaser of the product is permitted to copy this publication, for internal distribution only, from the media provided by Welch Allyn. No other use, reproduction, or distribution of this publication, or any part of it, is permitted without written permission from Welch Allyn.

Welch Allyn assumes no responsibility for any injury to anyone, or for any illegal or improper use of the product, that may result from failure to use this product in accordance with the instructions, cautions, warnings, or statement of intended use published in this manual.

Software in this product is Copyright 2019 Welch Allyn or its vendors. All rights are reserved. The software is protected by United States of America copyright laws and international treaty provisions applicable worldwide. Under such laws, the licensee is entitled to use the copy of the software incorporated with this instrument as intended in the operation of the product in which it is embedded. The software may not be copied, decompiled, reverse-engineered, disassembled, or otherwise reduced to human-perceivable form. This is not a sale of the software or any copy of the software; all right, title, and ownership of the software remain with Welch Allyn or its vendors.

For information about any Welch Allyn product, contact your local Welch Allyn representative: <http://www.welchallyn.com/en/other/contact-us.html>.

729098

This manual applies to REF 901110 NONMASTED SCALES.

DIR 80024375 Ver. A Revision date: 2019-01

Welch Allyn, Inc.
4341 State Street Road
Skaneateles Falls, NY 13153 USA

www.welchallyn.com

WelchAllyn[®]
Advancing Frontline Care[™]

Contents

Introduction	1
Intended use	1
Indications for use	1
Intended clinical care environments	1
Symbols	1
About warnings and cautions	3
Controls and indicators	5
Assembly and setup	7
Assembly instructions	7
Scale operation	9
Power options	9
Audible beeps	10
Weighing procedure	10
Diaper weighing procedure	11
Send measurements to a computer	12
Send measurements to a Welch Allyn device	13
Scale transportation	14
Storage, cleaning, and maintenance	15
Storage	15
Cleaning and disinfecting	15
Battery replacement	16
Maintenance	16
General compliance and standards	18
Troubleshooting	18
Appendices	21
Custom setup	21
Specifications	22
Product configurations	24
Approved accessories	24
Warranty	25
EMC compliance	25

Introduction

Intended use

The Welch Allyn Pediatric/Infant Scale is intended to be used by clinicians for weighing patients from 0.35 oz to 44 lbs (10 g to 20 kg) and measuring patients up to 32 inches (81 cm), depending on the size of the scale cradle.

Indications for use

The Welch Allyn Pediatric/Infant Scale is used by clinicians to weigh and measure the length of neonate and pediatric patients.

Pediatric scales can make contact with a patient's head, neck, back, arms, legs, and side. Contact duration is intended to be limited to less than 30 seconds.

Intended clinical care environments

The Welch Allyn Pediatric/Infant Scales are intended to be used in the following clinical care environments:

- Hospitals
- Ambulatory care centers
- Physicians' offices
- Other professional medical facilities

Symbols

Documentation symbols

WARNING The warning statements in this manual identify conditions or practices that could lead to illness, injury, or death.

Note

Warning symbols will appear with a grey background in a black and white document.

Caution The caution statements in this manual identify conditions or practices that could result in damage to the equipment or other property, or loss of data.

Mandatory - Consult Directions for Use

Shipping, storing, and environment symbols

Temperature limitation

Keep dry

Separate collection of Electrical and Electronic Equipment. Do not dispose as unsorted municipal waste.

Fragile

Humidity limitation

Atmospheric pressure limitation

This way up

Stacking limit by number

Recyclable

Keep away from sunlight

Atmospheric pressure limitation

Date of manufacture

Miscellaneous symbols

Manufacturer

Reorder Number

Product Identifier

For indoor use only

Serial Number

Prescription only or "For Use by or on the order of a licensed medical professional"

Battery

Direct current (DC)

	USB		On/Off Pushbutton
	Class II equipment		Type B applied part Note The entire scale is considered an applied part.
	Do not re-use		Mass
	Maximum safe working load limits		Rated power input, DC
	With respect to electrical shock, fire, and mechanical hazards only In accordance with: ANSI/AAMI ES60601-1:2005/(R)2012 and A1:2012, C1:2009/(R)2012 AND A2:2010/(R)2012, CAN/CSA C22.2 No. 60601-1:14 IEC 60601-1 Ed. 3.1		Lot code
	Global Trade Item Number		

About warnings and cautions

Caution statements can appear on the Welch Allyn Pediatric/Infant Scale, the packaging, the shipping container, or in this *Directions for use*.

The Welch Allyn Pediatric/Infant Scale is safe for patients and trained clinicians when the scale is used in accordance with the instructions and caution statements presented in this *Directions for use*.

Before using the device, you must familiarize yourself with all cautions, with the steps to power up the device, and with the sections of this *Directions for use* that pertain to your use of the device. In addition to reviewing the general cautions presented in the next section, you must also review the more specific cautions that appear throughout the manual in conjunction with setup/startup, operation, and maintenance tasks. No additional training is required.

- Failure to understand and observe any warning statement in this manual could lead to patient injury or illness.
- Failure to understand and observe any caution statement in this manual could lead to damage to the equipment or other property, or loss of patient data.

General warnings and cautions

WARNING Patient injury risk. Never leave the infant unattended on the scale. Only remove your hands from the infant for a brief amount of time to allow for an accurate weight measurement.

WARNING Patient or operator injury risk. Make sure the scale is positioned on a stable surface and in a way to avoid damage or bumping hazards.

WARNING Patient injury risk. Make sure the scale is set in the proper measurement units for your facility. Make sure you record the measurement units as they are displayed.

WARNING Electric shock hazard. Use only a Welch-Allyn approved power supply. The use of an unapproved power supply could increase chassis or patient leakage currents.

WARNING Electric shock hazard. Use only a Welch-Allyn approved power supply. The use of an unapproved power supply could cause electric harm and shock to you or the patient.

WARNING Electric shock hazard. All signal input and output (I/O) connectors are intended for connection of only devices complying with IEC 60601-1, or other IEC standards (for example, IEC 60950), as applicable to the scale. Connecting additional devices to the scale could increase chassis or patient leakage currents. To maintain operator and patient safety, consider the requirements of IEC 60601-1-1. Measure the leakage currents to confirm that no electric shock hazard exists.

WARNING Electric shock hazard. Make sure that the USB (SIP/SOP) port and the patient are never touched or come in contact at the same time.

WARNING No modification of this equipment is allowed.


CAUTION Make sure that you routinely perform general maintenance and equipment safety checks on your scale. Remove the scale from service when you notice damaged to the power cable.

CAUTION Do not use this scale to transport patients or items.


Controls and indicators

Scale

-
- 1 Weighing cradle with length gauge
 - 2 Rear bumper/cord wrap
 - 3 Carrying handles
 - 4 Front panel with push buttons and display
 - NA Battery cover (on top of scale)
 - NA DC power input (on rear of scale)
 - NA USB port (on rear of scale)
-

Front panel buttons

Button	Function
ZERO	Press this button to reset the scale value to zero.
Green LED indicators	This area displays the patient measurement, as well as indicates the measurement units, AC power, low battery, prior weight, and data transmission.
Power button 	Press this button to turn the scale on. Press and hold this button to turn the scale off.
SEND DATA	Pressing this button when the scale is connected to a computer sends the current measurement to the computer, a Welch Allyn Connex Vital Signs® Monitor, or a Welch Allyn Connex® Spot Monitor.
LB / KG	Press this button to change the scale between kilograms, pounds, ounces, and pounds and ounces.
REWEIGH	Press this button to perform a new weighing cycle.
RECALL	Press this button to display the weight stored in the scale's memory.

Assembly and setup

The scale is fully assembled at the factory and tested. In some instances, the scales are shipped assembled. However, due to size and shipping restrictions, most scales are partially disassembled. If assembly is required, it can be done with common hand tools. When special tools are required, they are supplied with the scale. Ensure that all bolts and fasteners are securely tightened.

Use the following instructions to assemble and setup your scale for use.

Assembly instructions

The scale is fully assembled at the factory and tested. The scale is shipped with the weighing cradle removed.

1. Carefully remove the acrylic weighing cradle from the shipping carton. If you ordered it, the oversized cradle ships in a separate carton.
2. Remove the scale from the shipping carton and place it on a stable surface. Make sure the front panel faces you.
3. Remove the two screws on the battery cover on the top of the scale, and remove the cover.
4. Insert the six D-size batteries that were included with your scale. Make sure to follow the polarity directions indicated on the battery holder.
5. Replace the battery cover, and re-install the two screws.
6. Position the cradle over the scale so the four metal mounting posts line up with the holes in the top of the scale.

Note the length gauge and position the cradle so the length gauge begins on the left.

7. Remove the protective cover from the cradle.
8. Apply gentle downward pressure to seat all four posts so the cradle rests firmly on the scale.
9. If you purchased the optional 4802D Pediatric Scale Cart, follow the instructions included with the cart.

Scale operation

Power options

The scale can be used with battery power or with AC power.

Press the power button to turn on the scale. When connected to AC power, the green AC power indicator illuminates.

Battery power

CAUTION Use only size D disposable alkaline batteries. The use of any other battery will void the warranty.

Press the power button to turn on the scale.

The scale automatically turns off after 90 seconds of non-use when operating on battery power.

The **LOW BATTERY** indicator illuminates when the battery power is low.

bAttrY displays when the batteries are unable to properly power the scale.

Replace the batteries when either of these messages are displayed.

AC power

When plugged in, the scale continuously runs on AC power automatically, depending on your settings.

To disconnect the scale from AC power, either remove the power adaptor from the wall or the back of the scale. Make sure to position the scale so that it can be easily disconnected from AC power.

The scale switches to battery power when AC power is not present.

Audible beeps

When the beeper option is on, the following beeps occur when using the scale.

Note The beeper option is set to OFF by default. These beeps do not occur unless beeper option has been enabled in the setup menu. See the Custom setup section to enable the beeper option.

Sound	Meaning
One beep	Acknowledgment of a button press.
Two beeps	A zero weight is obtained.
Three beeps	Weight reading is obtained and displayed.
Four long beeps	The battery is too weak to operate the scale. The scale has entered custom setup mode.
Long/short beeps	There is a problem with the scale. See the troubleshooting section.
Continuous long beep	There is a problem with the scale. See the troubleshooting section.

Weighing procedure

CAUTION Replace the cradle after any sign of damage, including small cracks or crazing. Never place a patient on a cradle exhibiting signs of damage.

Make sure to select a consistent time of day to obtain a weight measurement.

Position the scale in a convenient and stable location.

WARNING Patient injury risk. Make sure the scale is properly cleaned between patients.

1. Press the power button to turn on the scale.

The scale goes through the startup sequence.

2. Prepare the infant for weighing.
3. Place a scale liner on the scale cradle.

WARNING Patient injury risk. Do not reuse the scale liners. They are single-patient use only.

4. Press **ZERO** to zero the scale.
5. When the scale display reads **0**, place the infant on the scale so that the weight is distributed evenly near the center of the scale.

WARNING Patient injury risk. Never leave the infant unattended on the scale. Only remove your hands from the infant for a brief amount of time to allow for an accurate weight measurement.

6. Stabilize the infant as necessary with your hand. Remove your hand from the infant, but always keep your hands close in case re-stabilization is required.

Note Do not touch the infant while the weight measurement is in process.

7. The scale displays a series of moving dashes as it waits for stabilization and calculates the weight.

Once the measurement is obtained, the scale beeps three times and displays the weight on the front panel.

Note The beeper option is set to OFF by default. These beeps do not occur unless beeper option has been enabled in the setup menu. See the Custom setup section to enable the beeper option.

WARNING Patient injury risk. Make sure the scale is set in the proper measurement units for your facility. Make sure you record the measurement units as they are displayed.

8. If necessary, press **REWEIGH** to perform another weight measurement.
9. Return the infant to his or her bassinet or incubator.
10. If necessary, press **LB/KG** to change the unit of measure on the display.

Note If you purchased the kilogram-only option, the scale displays **Lb OFF** when **LB/KG** is pressed.

11. If necessary, press **RECALL** to display the previously stored scale measurement.

The **PRIOR WEIGHT** LED flashes.

Note A measurement is stored in the scale until a new weight is taken and zeroes out.

A series of moving dashes indicates that a tared weight is stored in the scale memory, or that there is a negative weight on the scale. Press **ZERO** to clear the tare value and reset the scale.

If **O-load** is displayed, the capacity of the scale has been exceeded.

Diaper weighing procedure

Use the following procedure to determine urine output using diaper weight. When displaying weights between 10 grams to 1 kilogram, the scale automatically converts to a 1-gram (0.001 kilogram) resolution.

Make sure that the scale is located in still air, and away from breezes caused by fans or HVAC systems.

1. Press **ON** to turn on the scale.
2. (Optional) Press **ZERO** to zero the scale.
3. Place a new, dry diaper of the identical size and brand as the diaper to be weighed on the weighing platform, and make sure that the scale is in kilograms mode.

Three beeps sound and the weight is displayed.

Note The beeper option is set to OFF by default. These beeps do not occur unless beeper option has been enabled in the setup menu. See the Custom setup section to enable the beeper option.

4. After the weight is displayed, press **ZERO**.

0000 is displayed.

5. Replace the dry diaper with the wet diaper on the center of the cradle.

The scale displays a series of moving dashes until the wet diaper is detected. Once detected, three beeps sound and the weight displays.

Note The beeper option is set to OFF by default. These beeps do not occur unless beeper option has been enabled in the setup menu. See the Custom setup section to enable the beeper option.

6. The displayed weight is the difference between the wet and dry diaper. Record the measurement displayed on the scale.

Note One gram (0.001 kg) is equal to 1 cc of liquid. (1.0 specific gravity).

To ensure the most accurate measurements, do the following:

- Perform the procedure each time a urine output measurement is needed.
- Do not rely on subtracting a known weight for a dry diaper from the displayed weight of a wet diaper.
- Remove any fecal matter from the wet diaper before starting the weighing procedure. The change in the dry and wet diaper should only reflect urine content.

A series of moving dashes indicates that a tared weight is stored in the scale memory, or that there is a negative weight on the scale. Press **ZERO** to clear the tare value and reset the scale.

Note When using battery power, the scale does not save tare values after powering down. See the Custom setup section to change the number of seconds before the scale turns off when operating on battery power.

Send measurements to a computer

The scale can connect to a computer using the Microsoft® Windows® 10 operating system.

Detailed instructions for this feature are available from Welch Allyn. Go to <http://www.welchallyn.com/en/other/contact-us.html> to find your local representative.

CAUTION Any computer connected to the scale must be running on battery power, a 60601-1-compliant power supply, or a 60601-1-compliant isolation transformer.

1. Connect the USB Type B male side of the cable to the USB port on the rear of the scale.
2. Connect the USB Type A male side of the cable into a USB port on your computer. Windows automatically installs the necessary driver.

3. Press **SEND DATA** on the scale to send the current measurement to the computer.

Send measurements to a Welch Allyn device

The scale can connect to a Welch Allyn Connex[®] Vital Signs Monitor or a Welch Allyn Connex[®] Spot Monitor.

CAUTION Any device connected to the scale must be running on battery power, a 60601-1-compliant power supply, or a 60601-1-compliant isolation transformer.

1. Connect the USB Type B male side of the cable to the USB port on the rear of the scale.
2. Connect the USB Type A male side of the cable into a USB port on the rear of your device.
3. Press **SEND DATA** on the scale to send the current measurement to the device.

Scale transportation

Manual transport

Note If your scale has an oversize cradle, make sure to remove and safely handle the cradle before moving the scale.

1. Unplug the power adapter cable from the scale.
2. Use only the handles to lift and move the scale.
3. At the new scale location, do the following:
 - a. Make sure the scale is positioned on a flat surface with the rear bumper facing the wall.
 - b. Make sure the scale is positioned so that all four feet are on a stable surface.
 - c. If you have an oversize cradle, make sure to place the cradle on the scale properly.
 - d. Make sure to plug the power adapter cable into the wall and into the back of the scale.
 - e. Make sure the power adapter cable does not pose a tripping hazard and is easily accessible.

Cart transport

The cart is a purchasable option.

1. Make sure the cart wheels are unlocked.
2. Unplug the power adapter cable from the rear of the scale.
3. Use only the cart handle to move the scale.
4. At the new scale location, do the following:
 - a. Make sure the scale is positioned so that all four feet are on top of the cart.
 - b. Make sure to lock the wheels on the cart.
 - c. Make sure to plug the power adapter cable into the wall and into the back of the scale.
 - d. Make sure the power adapter cable does not pose a tripping hazard and is easily accessible.

Storage, cleaning, and maintenance

Storage

The scale should be stored on a stable surface. If not, the scale should be stored in a convenient storage facility or closet. At no time should the scale be lifted by grasping the cradle. During storage or shipping, disconnect the cradle from the scale to prevent damage to the scale's electronic transducers.

Cleaning and disinfecting

Scale

Note Any visible soil must be removed before cleaning or disinfecting.

1. Use one of the following products when cleaning or disinfecting your scale:
 - PDI Super Sani-Cloth® (disinfectant) and Sani-Cloth® Plus (disinfectant)
 - Metrex CaviWipes™ (disinfectant)
 - A clean cloth moistened with 70 percent isopropyl alcohol
 - Clorox Healthcare® Bleach Germicidal Wipes (disinfectant) or Sani-Cloth® Bleach Wipes (disinfectant)
2. For optimal results, follow the manufacturer's instructions and your facility's procedure to clean the scale. Make sure to properly dispose of any used cleaning supplies.

CAUTION The scale will tarnish with prolonged contact with bleach. Use only the approved chemicals indicated above.

Cradle

Note Any visible soil must be removed before cleaning or disinfecting.

1. Use one of the following products when cleaning or disinfecting the cradle:
 - PDI Sani-Cloth® AF3 (disinfectant)
 - Metrex CaviWipes™ AF (disinfectant)
 - A clean cloth moistened with mild soap and water
 - PDI Sani-Cloth® Bleach Germicidal Disposable Wipe (disinfectant)

2. For optimal results, follow the manufacturer's instructions and your facility's procedure to clean the scale. Make sure to properly dispose of any used cleaning supplies.
3. (Optional) Perform the following steps to remove deposits of dissolved solids after the disinfecting liquid evaporates:
 - a. Wipe the cradle with a soft, clean cloth moistened with clean water.
 - b. Dry the cradle with a soft, dry cloth.

CAUTION Do not clean the cradle with alcohol or other drying agents.

Battery replacement


WARNING Electronic shock hazard. Do not replace the battery under the following conditions:

- in the vicinity of the patient
- with wet hands
- in the presence of flammable anesthetics

CAUTION Use only size D disposable alkaline batteries. The use of any other battery will void the warranty.

1. Make sure the scale is powered down.
2. Remove the cradle from the scale.
3. Remove the battery cover by removing the two screws with a Phillips screwdriver.
4. Remove the old batteries and install six new size D batteries in the battery holder. Make sure to follow the polarity instructions.

5. Re-attach the top cover of the enclosure by installing the two screws.
6. Re-install the cradle by gently pressing down.

Maintenance

Perform the following preventive maintenance quarterly to keep your scale in working order, or as required based on usage or per your facility maintenance schedule, whichever comes first.

1. Have your service department check the calibration annually or as required. Calibration instructions are provided in the scale service manual.
2. Inspect the cradle for cracks or loose mounting hardware. Replace or repair as necessary.
3. Inspect the scale's enclosure for damage or loose or missing hardware. Replace or repair as necessary.
4. Inspect the power adapter cable for abrasions or other signs of wear.
5. Do not expose the scale to excessive water or moisture.

6. Do not store the scale where heavy objects can be placed on it.
7. Replace the batteries annually or as required.
8. When storing the scale, remove the batteries from the battery area. Batteries can corrode over a period of time. Make sure to check the batteries before putting the scale back into use.
9. Do not service or perform any maintenance while the scale is in use with a patient.

General compliance and standards

**Directive 2002/96/EC-WEEE:
Disposal of noncontaminated electrical and electronic equipment**

This product and its accessories must be disposed of according to local laws and regulations. Do not dispose of this product as unsorted municipal waste. Prepare this product for reuse or separate collection as specified by Directive 2002/96/EC of the European Parliament and the Council of the European Union on Waste Electronic and Electrical Equipment (WEEE). If this product is contaminated, this directive does not apply.

For more specific disposal or compliance information, see www.welchallyn.com/weee, or contact Welch Allyn Customer Service at +44 207 365 6780.

Troubleshooting

This section presents tables of information to help you troubleshoot issues with the scale.

To use these tables, locate the specific problem with your scale in the left column of the table. The remainder of the row explains possible causes and suggests actions that can resolve the issue.

Issue	Possible cause	Suggested action
The display does not illuminate	No power is supplied to the scale	Make sure that the scale is plugged in. Make sure that the batteries are properly installed and are in good condition.
	The cable between the display board and the instrument board is disconnected	Contact your service department.
	The membrane keyboard and the display board are disconnected	Contact your service department.
The weight reading takes too long to display	The cradle is in motion	Wait for the cradle to settle.
The LOW BATTERY indicator is illuminated	The scale batteries are low on power	Replace the scale batteries.
The scale displays bAttrY	The batteries do not have enough power to operate the scale	Replace the scale batteries.
The scale displays CABLE	A signal outside the range of the internal A/D converter has been applied	Contact your service department.
The scale displays O-LOAD	The object on the scale is larger than the maximum value allowed	If the weight value is within the specified range, contact your service department.

Issue	Possible cause	Suggested action
The reading does not change when weight is applied to the scale	The load cell transducer, connector, or cable may be defective	Contact your service department.
The scale buttons are not responding	The buttons are damaged	Check for visible signs of damage. Contact your service department.
The error message ERR FD is displayed	The scale's factory database has failed	Contact your service department.
The error message ERR DB is displayed	The scale's application database has failed	Contact your service department.
The error message ERRoR is displayed	The scale's flash memory has failed	Contact your service department.
The error message ERR AD is displayed	The measurement board cannot be read	Contact your service department.
The error message ERR AV is displayed	The measurement board is the wrong version	Contact your service department.
The error message ERR CL is displayed	Measurement board error, the weight display is blocked	Contact your service department.

Appendices

Custom setup

You can customize the scale to best suit your needs. Configurations that can be selected include the automatic shut-off time, weighing units, beeper, and more.

Enter custom setup mode

1. Make sure that the scale is powered down.
2. Press and hold **ZERO** and **REWEIGH** while pressing the power button.
0000 is displayed.
3. Enter the PIN **9821** using the following keys:
 - Press **REWEIGH** to adjust the selected value up or press **RECALL** to adjust the selected value down.
 - Press **LB/KG** to confirm the entry and move to the next value.
4. Once the correct PIN is entered, press **LB/KG** to enter the setup mode
The scale displays **Set-UP** and produce a series of four long beeps.

Set the options

Enter the custom setup as described.

Press **LB/KG** to scroll through the custom setup options.

The options displayed with each press are as follows:

Press **REWEIGH** to adjust these options up or press **RECALL** to adjust these options down. Press **LB/KG** to accept the current value and move on.

Note Options indicated with an asterisk (*) require an additional press of **LB/KG** to change the value.

Option displayed	Feature
SOft	This option displays the software version of your scale.
dAtE	This option displays the release date of the software. The format is MM.DD.YY.
SCALE	This option displays the model number of the scale.

AutOFF*	This displays the number of seconds before the scale turns off when operating on battery power. The Cont value prevents the scale from turning off automatically.
ACCont*	This option causes the scale to remain on when plugged in to AC power.
rES	Pressing REWEIGH or RECALL switches between the following options: 0.001 kg 0.002 kg 0.005 kg 0.010 kg
UnitS*	This option allows you to change the displayed weight unit. Note Do not change the scale units if you have purchased the kilogram-only option. KILOS* This option allows the weight to display in kilograms. POUNDS and OUNCES* This option allows the weight to display in pounds and ounces. POUNDS* This option allows the weight to display in pounds. OUNCES* This option allows the weight to display in ounces.
bEEPEr*	This option determines the audible signal that occurs when a front panel button is pressed, a weight is obtained, and other actions occur. By default, this option is set to OFF.
rEcALL	This option allows you to turn on the scale in Recall mode. If this option is set to OFF, the RECALL button does not operate.
USb*	This option allows you to send weight values to a computer.
Prtunt*	This option determines whether pounds or kilograms are used in the send data function when the scale is connected to a computer or a Welch Allyn device. Press REWEIGH or RECALL to switch between the following options: If the pounds and kilograms indicators are illuminated on the control panel, the data sent is determined by whichever unit is selected on the control panel. If the kilogram indicator is illuminated, it will only send data in kilograms, regardless of the unit selected on the control panel. If the pounds indicator is illuminated, it will only send data in pounds, regardless of the unit selected on the control panel.
PrtOPt*	Press REWEIGH or RECALL to switch between the following options: PnlPrt : The current weight is sent to a computer or a Welch Allyn device when SEND DATA is pressed on the control panel. AutPrt : The current weight is automatically sent when a weight reading occurs.
SEt-UP	This option is displayed when you have cycled through all the options. Press and hold the power button to power down the scale.

Specifications

This device is powered by an external power supply. When used together, the device and the power supply are considered ME equipment.

Item	Specification
Scale	4802D
Accuracy	5 g for weights between 10 g to 1 kg 0.5% for weights between 1 kg to 20 kg
Resolution (factory default)	$\pm .001$ kg from 10 to 1000 g and $\pm .005$ kg for 1000 g and above
Weighing range	0.35 oz to 44 lbs (10 g to 20 kg)
Environmental specifications	<p>Operation</p> <ul style="list-style-type: none"> • Temperature : 10 °C to 40 °C (50 °F to 104 °F) • Humidity: 10% RH to 90% RH, non-condensing. • Atmospheric pressure: 700 hPa to 1060 hPa <p>Shipping</p> <ul style="list-style-type: none"> • Temperature: -20 °C to 50 °C (-4 °F to 122 °F) • Humidity: 15% RH to 95% RH • Atmospheric pressure: 700 hPa to 1060 hPa <p>Storage</p> <ul style="list-style-type: none"> • Temperature -35 °C, humidity uncontrol • Temperature 30 °C, humidity 90% • Temperature 55 °C, humidity 30%
Protection rating	IPX0
Readout	Large, bright, light-emitting diode display. Selectable weight display in ounces, kilograms, pounds, or pounds and ounces. Indicator light to show selected units.
Power source	<p>Battery</p> <ul style="list-style-type: none"> • 6 D alkaline batteries (1.5V x 6) <p>Power adapter:</p> <ul style="list-style-type: none"> • Medical grade UL listed external power supply • Nominal input voltage: 100-240V AC $\pm 10\%$ • Nominal input frequency: 50-60 Hz • Nominal input current: 0.4-0.2Arms @ max load • Standby power consumption at U_{in}: 230V AC : $\leq 0.1W$ • Nominal output volt age: U_{out} : 12V DC +5% / -5% U_{Br} : $\leq 150mVpp$ • Nominal output current: I_{out} : 1500mA
Dimensions	<p>Platform with carrying handles: 27.5 in x 16 in. x 7 in.</p> <p>Cradle:</p> <ul style="list-style-type: none"> • Standard – 25 in. x 14 in. • Oversized – 32 in. x 14 in. • Four Sided – 25 in. x 14 in.
Scale weight	27 lbs (12.2 kg) with the oversize cradle

Technical description

The Welch Allyn Pediatric/Infant Scale uses load cell transducers to convert analog weight applied to the weighing platform into a digital, electrical signal for display on the readout.

Welch Allyn will provide, on request, circuit diagrams, component part lists, descriptions, calibration instructions, and other information to authorized service personnel

Product configurations

The Welch Allyn Pediatric/Infant Scale is available in the following configurations:

Model 4802D, Pediatric/Infant Scale

Model	Description
4802D-AK-FB	4802D Pediatric Scale / Infant Scale with four-sided cradle, metric units only (kg)
4802D-AK-VB	4802D Pediatric Scale / Infant Scale with oversized cradle, metric units only (kg)
4802D-AK-XB	4802D Pediatric Scale / Infant Scale with standard cradle, metric units only (kg)
4802D-AX-FB	4802D Pediatric Scale / Infant Scale with four-sided cradle, imperial and metric units (lb/kg)
4802D-AX-VB	4802D Pediatric Scale / Infant Scale with oversized cradle, imperial and metric units (lb/kg)
4802D-AX-XB	4802D Pediatric Scale / Infant Scale with standard cradle, imperial and metric units (lb/kg)

Approved accessories

Part number	Description
412490	4802D Pediatric Scale Cart
726551	Power supply
48200W	Standard cradle assembly
48220W	Four-sided cradle assembly
48222	Oversize cradle assembly

Additional items

Part number	Description
48217	Scale liners (500 per carton, 100 per bag)

Warranty

Welch Allyn will warranty the weight scale to be free of defects in material and workmanship and to perform in accordance with manufacturer specifications for the period of one year from the date of retail purchase.

The warranty period shall start on the date of purchase. The date of purchase is: 1) the invoiced ship date if the device was purchased directly from Welch Allyn, 2) the date specified during product registration, 3) the date of purchase of the product from a Welch Allyn authorized distributor as documented from a receipt from said distributor.

This warranty does NOT cover damages caused by misuse or abuse, including but not limited to:

- Failure caused by unauthorized repairs or modifications
- Damage caused by shock or dropping during transportation
- Damage caused by improper use of the power supply
- Failure caused by improper operation not consistent with the instructions stated in this *Directions for use*

Should this device require maintenance (or replacement at our option) under warranty, contact your local Welch Allyn representative: <http://www.welchallyn.com/en/other/contact-us.html>

EMC compliance

Special precautions concerning electromagnetic compatibility (EMC) must be taken for all medical electrical equipment. This device complies with IEC 60601-1-2 Ed 4.0 / EN 60601-1-2:2015.

- All medical electrical equipment must be installed and put into service in accordance with the EMC information provided in this document, the *Welch Allyn® Pediatric/Infant Scale Directions for use*.
- Portable and mobile RF communications equipment can affect the behavior of medical electrical equipment.

The device complies with all applicable and required standards for electromagnetic interference.

- It does not normally affect nearby equipment and devices.
- It is not normally affected by nearby equipment and devices.
- It is not safe to operate the scale in the presence of high-frequency surgical equipment.
- However, it is good practice to avoid using the device in extremely close proximity to other equipment.

Emissions and immunity information

Electromagnetic emissions

The Welch Allyn Pediatric/Infant Scale is intended for use in the electromagnetic environment specified below. The customer or user of the Welch Allyn Pediatric/Infant Scale should assure that it is used in such an environment.

Electromagnetic emissions

Emissions test	Compliance	Electromagnetic environment - guidance
Radiated emissions CISPR 11	Group 1 /Class B	The Welch Allyn Pediatric/Infant Scale uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
Conducted emissions CISPR 11	Class B	The Welch Allyn Pediatric/Infant Scale is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Electromagnetic immunity

The Welch Allyn Pediatric/Infant Scale is intended for use in the electromagnetic environment specified below. The customer or the user of the Welch Allyn Pediatric/Infant Scale should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode	±1 kV differential mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<0 % U_T (>100 % dip in U_T) for 1 cycle @ 0° 70 % U_T (30 % dip in U_T) for 25/30 cycles @ 0° <0 % U_T (>100 % drop in U_T) for 0.5 cycles @ 0°, 45°, 90°, 135°, 180°, 225°, 270°, & 315° 0% for 250/300 cycles @ 0°	<0 % U_T (>100 % dip in U_T) for 1 cycle @ 0° 70 % U_T (30 % dip in U_T) for 25/30 cycles @ 0° <0 % U_T (>100 % drop in U_T) for 0.5 cycles @ 0°, 45°, 90°, 135°, 180°, 225°, 270°, & 315° 0% for 250/300 cycles @ 0°	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Welch Allyn Pediatric/Infant Scale requires continued operation during power mains interruptions, it is recommended that the Welch Allyn Pediatric/Infant Scale be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Electromagnetic immunity

Note: U_T is the a.c. mains voltage prior to application of the test level.

Electromagnetic immunity

The Welch Allyn Pediatric/Infant Scale is intended for use in the electromagnetic environment specified below. The customer or the user of the Welch Allyn Pediatric/Infant Scale should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
			Portable and mobile RF communications equipment should be used no closer to any part of the Welch Allyn Pediatric/Infant Scale, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
			Recommended separation distance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 V	$d = (1.17) \sqrt{P}$
	6Vrms ISM Band frequency (MHz) : 6.675 – 6.795, 13.553 – 13.567, 26.957 – 27.284 & 40.66 – 40.70	6V	
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz	3 V/m	$d = (1.17) \sqrt{P}$ 80 to 800 MHz

$$d = (2.33) \sqrt{P} \text{ 800 MHz to 2.7 GHz}$$

where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range^b. Interference may occur in the vicinity of equipment marked with the following symbol:

Table-9 of IEC60601-1-2, 4th Ed.

27 V/m	27 V/m
380 MHz to 390 MHz	380 MHz to 390 MHz
28 V/m	28 V/m

Electromagnetic immunity

430 MHz to 470 MHz	430 MHz to 470 MHz
800 MHz to 960 MHz	800 MHz to 960 MHz
1700 MHz to 1990 MHz	1700 MHz to 1990 MHz
2400 MHz to 2570 MHz	2400 MHz to 2570 MHz
9 V/m	9 V/m
704 MHz to 787 MHz	704 MHz to 787 MHz
5100 MHz to 5800 MHz	5100 MHz to 5800 MHz

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Welch Allyn Pediatric/Infant Scale is used exceeds the applicable RF compliance level above, the Welch Allyn Pediatric/Infant Scale should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Welch Allyn Pediatric/Infant Scale.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the Welch Allyn Pediatric/Infant Scale

The Welch Allyn Pediatric/Infant Scale is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the Welch Allyn Pediatric/Infant Scale can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Welch Allyn Pediatric/Infant Scale as recommended below, according to the maximum output power of the communications equipment.

Separation distance according to frequency of transmitter (m)

Rated max. output power of transmitter (W)	150 kHz to 80 MHz $d = (1.17) \sqrt{P}$	80 MHz to 800 MHz $d = (1.17) \sqrt{P}$	800 MHz to 2.7 GHz $d = (2.33) \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.20	1.20	2.30
10	3.70	3.70	7.40
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Recommended separation distances between portable and mobile RF communications equipment and the Welch Allyn Pediatric/Infant Scale

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.


WelchAllyn[®]
Advancing Frontline Care[™]

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