



Hillrom™

CERTIFIED PREVENTIVE MAINTENANCE



Beds and Surfaces

REDUCE SAFETY RISKS AND STAY COMPLIANT WITH:

✓	Documentation <ul style="list-style-type: none">■ Accreditation Canada Requirements
✓	Electrical Safety Certification <ul style="list-style-type: none">■ Measure leakage and ground resistance■ Inspect power cord and line filter for damage
✓	Functional Controls <ul style="list-style-type: none">■ Ensure caregiver and patient controls are operating correctly
✓	Inspecting Internal Components <ul style="list-style-type: none">■ Air blowers, mattress bladders and hoses■ Check mattress interior for moisture or stains
✓	Scale Calibration
✓	Inspecting Safety Features <ul style="list-style-type: none">■ Siderails■ CPR feature■ Alerts, SideCom® Communication System, Bed Exit, SafeView® + System and wireless connectivity
✓	Inspecting/Adjusting Movement Limits <ul style="list-style-type: none">■ Trendelenburg, Auto Contour™ feature, chair, hi-low, head articulation, foot retraction/articulation, and SlideGuard® feature
✓	Inspecting Bed Movement Functions <ul style="list-style-type: none">■ Brakes, steering and 5th wheel■ Powered drive system
✓	Battery Backup <ul style="list-style-type: none">■ Check that battery activates when bed is unplugged and has sufficient power■ Confirm caregiver controls function with battery backup
✓	Improving Appearance <ul style="list-style-type: none">■ Replace labels and plastic parts

hillrom.ca

* The tests listed are typical and customary. All tests may not be applicable to all equipment.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2022 Hill-Rom Services, Inc. ALL RIGHTS RESERVED. 168813 rev 4 23-MAY-2022 ENG – CA

Hillrom is a part of Baxter

Baxter