## Welch Allyn Connex® ProBP 3400 Converting Handheld to Other Configurations

## To convert to Wall Mount Configuration, customer will need:

PN: 4601-61 Desc: WALL MOUNT WITH BASKET FOR PROBP 3400

PN: 3400-31 Desc: DOUBLE TUBE BLOOD PRESSURE HOSE (10FT) - OPTIONAL (to extend

the BP hose from 5' as supplied with the handheld configuration to 10' for extra

reach)



## To convert to Desk Mount Configuration, customer will need:

PN: 4602-61 Desc: DESK MOUNT W/BASKET FOR PROBP 3400



## To convert to Mobile Stand Configuration, customer will need:

PN: 4600-61 Desc: MOBILE STAND W/MOUNTING KIT PROBP 3400

PN: PWCD-5WT-B

Desc: PROBP MOBILE STAND CORD ASSY - OPTIONAL (4600-61 comes with a Power Supply Mounting Kit. The PWCD-5W-B includes a new transformer that mounts in

the Power Supply Mount, a short 18" USB cable that connects the output of the

transformer to the ProBP, and a long 8' power cord that plugs the transformer into a wall outlet. The power supply and long USB cable that comes with the Handheld configuration can be used with the mobile stand, but the PWCD-5WT-B provides a cleaner installation.)



