PROGRESSA® SMART BED SYSTEM

Empower patients and caregivers with in-bed therapies to help get your most critical patients better sooner.



//////

STANDARD FEATURES TO HELP OPTIMIZE PATIENT SAFETY AND IMPROVE CAREGIVER SATISFACTION

SlideGuard™ Frame Technology
Patient Controls Backlighting
FlexAfoot™ Retractable Foot Control
Advanced MicroClimate® Technology
30°/45° Head of Bed Alarm
In-Bed Scale
Line Managers
Three Mode Bed Exit Alarm
Caregiver Pendant
Caregiver Foot Controls
Night Light
Integrated Oxygen Tank Holder
Accessory Outlet
Hands Free Emergency CPR & Trendelenburg
NaviCare Ready - Local Alarm On/Off
Point of Care™ Siderails Control
Both Siderails with Point-of-Care Touch Screen
Smart Bed Ready (Sidecom Required)
Obstacle Detect™ System at Both Sides of the Bed
Drainage Bag Holders
Four IV Sockets

OPTIONAL FEATURES TO SUPPORT THERAPY-FOCUSED TECHNOLOGIES

Prevention S	Surface
Therapy Sur	face
Pulmonary S	Surface
Continuous	Lateral Rotation Therapy (CLRT)
Percussion &	& Vibration Therapy (P&V)
Nurse Call, L	Jniv TV & Lighting
Siderail Nurs	se Call
Bed Exit Alar	rm to Nurse Call
5th Wheel S	teering
IntelliDrive®	Powered Transport Steering
Transport Sh	nelf
Wireless Co	nnectivity
Permanent l	V Pole
Full Chair Eg	iress
StayInPlace™	MPatient Migration Management System
WatchCare"	^M Incontinence Management System

Progressa[®] Smart Bed System

TECHNICAL SPECIFICATIONS

Total Length Foot Section Extended Foot Section Retracted	98" (248.9 cm) 88" (223.5 cm)
Maximum Width Side Rails Down Side Rails Up	40.5" (103 cm) 40.5" (103 cm)
Maximum Side Rail Height Above Sleep Deck	17" (43.2 cm)
Highest Position Top of Seat Section to Floor	35.7" (90.7 cm)
Lowest Position Top of Seat Section to Floor	16.5" (41.9 cm)
Minimum Underbed Clearance	4.3" (10.9 cm) w/o IntelliDrive® Powered Transport
Standard Caster Size	6" (15.2 cm)
Articulations Head Section Thigh Section Foot Section	0° - 67° 77° with Chair Egress 0° - 30° 0° - 45°
Trendelenburg Position (maximum) Reverse Trendelenburg Position (maximum)	75° with Chair Egress 13° 18°

Safe Working Load	650 lb (295 kg)
Maximum Patient Weight	500 lb (227 kg)
C-Arm Clearance	7 - 15" (17.8-38.1 cm)
Cleanability/Infection Control	IPX4
Radiolucent Sleep Deck	Standard 17.7" x 23" (45 cm x 58.5 cm)
Voltages	100/110/115/120/ 127/220/230/240
Mattress Prevention Surface Mattress Width Mattress Length Mattress Height Mattress Weight Therapy Surface Mattress Width Mattress Height Mattress Weight Pulmonary Surface Mattress Width Mattress Length Mattress Length Mattress Height Mattress Height Mattress Weight	35" (88.9 cm) 84" (213 cm) 71" (18.1 cm) 31 lb (14.1 kg) 35.5" (90.2 cm) 84" (213 cm) 8" (20.32 cm) 45 lb (20.5 kg) 35.5" (90.2 cm) 84" (213 cm) 8" (20.32 cm) 48 lb (21.8 kg)

SURFACE & FRAME CONFIGURATION OPTIONS

PROGRESSA [®] PREVENTION	PROGRESSA® THERAPY	PROGRESSA [®] PULMONARY
Prevention Surface	Therapy Surface	Choose one: - Pulmonary Surface (CLRT only) - Pulmonary Surface (CLRT & P&V)
 Choose one: Dining Chair[®] Position FullChair[®] Position with Chair Egress FullChair[®] Position with Chair Egress and StayInPlace[™] Technology 	 Choose one: Dining Chair[®] Position FullChair[®] Position with Chair Egress FullChair[®] Position with Chair Egress and StayInPlace[™] Technology 	 Choose one: Dining Chair[®] Position FullChair[®] Position with Chair Egress FullChair[®] Position with Chair Egress and StayInPlace[™] Technology



hillrom.com

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.