

Tri-Flex™ II Bariatric Bed

Dependable Bariatric Care



Reliable and dependable choice for the care of bariatric patients.

Hospitals around the country face many challenges, including cutbacks, staffing shortages, and limited resources. Because of the unique medical demands, bariatric patients often require extra care, which can place even more strain on resources and staff.

The Tri-Flex II bariatric bed delivers:

- Dependable care
- Reliability
- Patient comfort

Hill-Rom can help **transform bariatric care.**

Enhancing outcomes for
patients and their caregivers:

Hill-Rom

Tri-Flex™ II Bariatric Bed

Dependable Bariatric Care

Technical Specifications

Bed

Low position—sleep deck to floor	14" (36 cm)
High position—sleep deck to floor	28" (71 cm)
Maximum overall width	40.75" (109 cm)
Maximum overall length	94" (939 cm)
Patient weight capacity	1,000 lbs (453.6 kgs)

Bed Angles

Head section	0-60°
Knee section	0-24°
Max Trend	12°

Sleep Surface

Sleep deck width	37"-48" (Rental) (93.98 cm-121.92 cm) 37"-54" (Capital) (93.98 cm-137.16 cm)
Sleep deck length	86" (137.16 cm)
Foam surface	86"L x 37"-48"W x 9.4"H (218.44 cm x 93.98-121.92 cm x 23.876 cm) [includes two 5.5" (13.97 cm interchangeable side bolsters)]
Low air loss surface (Synergy® Elite Surface)	86"L x 39" or 48"W x 10"H (218.44 cm x 99.06 cm or 121.92 cm x 25.4 cm)

Caster & Braking System

Caster diameter	6" (15.26 cm)
Four wheel braking	No
Steer and Lock	Yes

Scale System

Accuracy	1% of patient weight
Scale weight capacity	1,000 lbs (453.6 kgs)

Safety Listing and Power Requirements

Safety listing	UL 60601-1
Power requirements	120v AC , 60H

Standard Features

Expandable width frame and foam surface
Pendant
Battery backup

Optional Features

Full frame trapeze
In-bed scale (standard on rental)
Synergy® Air Elite surface
Impulse Drive System — Capital Purchase Only

Tri-Flex™ II is a trademark of Burke, Inc.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2015 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.
145598 rev 7 29-OCT-2015 ENG – US

For further information about this product or a service, please contact your local Hill-Rom representative or visit our webpage:

USA 800-445-3730
Canada 800-267-2337

www.hill-rom.com

Enhancing outcomes for
patients and their caregivers:

Hill-Rom