



## Solid State Headlights from Welch Allyn

FOR VETERINARIANS



High-quality, lightweight solid state headlights made affordable for every veterinary practice



### Solid State Procedure Headlight

- 50 lumens of bright, white, shadow-free light for true tissue color
- Cool beam, safe to touch
- At just six ounces, it's one of the most lightweight and comfortable headlights available
- Estimated life of up to 10,000 hours of continuous burn time
- Variable spot size for well-defined illumination
- Available with direct power source or the rechargeable "belt-clip" portable power source



### Solid State Portable Headlight

- 20 lumens of cool, white light for true tissue color
- Weighing just four ounces, it's one of the most lightweight headlights available
- Estimated life of up to 10,000 hours of continuous burn time
- Variable spot size for well-defined illumination
- Rechargeable battery housed in soft, terry cloth headband allows the user freedom of movement

**WelchAllyn**<sup>®</sup>

Advancing Frontline Care<sup>™</sup>



# Solid State Headlight Specifications

## Solid State Procedure Headlight

### Specifications

**Light Source:** 5-watt solid state lamp, color temperature 5500° Kelvin, 50 lumens, estimated life of up to 10,000 hours continuous burn

**Working Distance:** 16 inches

**Spot Size Range:** 1.5 to 6 inches at the working distance

**Weight:** 6 ounces

**Battery Run Time:** Approximately 1-hour continuous run time on a fully charged battery

**Battery Life:** Approximately 2.5 years (depending on usage)

**Warranty:** 1-year limited warranty

### Model Numbers

**49020** Solid State Procedure Headlight with Direct Power Supply

### Accessory Numbers

**49099** Carrying Case

**75260** Battery Pack with Charger

### Replacement Parts

**49000** Procedure Headlight only

**72250** Replacement Battery Pack

**74180** Charger

**75270** Direct Power Supply (U.S.)

## Solid State Portable Headlight

### Specifications

**Light Source:** 1-watt solid state lamp, color temperature 5500° Kelvin, 20 lumens, estimated life of up to 10,000 hours continuous burn

**Working Distance:** 16 inches

**Spot Size Range:** 1.5 to 6 inches at the working distance

**Weight:** 4 ounces

**Battery Run Time:** Approximately 90 minutes of continuous run time on a fully charged battery

**Battery Life:** Approximately 2 years (depending on usage)

**Warranty:** 1-year limited warranty

### Model Numbers

**46070V** Solid State Portable Headlight with Direct Power Supply/Charger

### Accessory Numbers

**05232-U** Soft Case

### Replacement Parts

**49642V** Soft Terry Cloth Headband—case of 10

**72270** Rechargeable Battery

**74170** Direct Power Supply (U.S.)



For international part numbers please contact Welch Allyn Customer Service or your local Welch Allyn distributor.

4341 State Street Road, PO Box 220, Skaneateles Falls, NY 13153-0220 USA  
 (p) 1.315.685.4560 or 1.800.535.6663 (f) 1.315.685.3361 [www.welchallyn.com/dvm](http://www.welchallyn.com/dvm)