


# Flexible Illumination Systems for Optimal Visualization



## Sigmoidoscope/Anoscope Lighting Systems

Powerful illumination systems.

- > Halogen light provides consistent, long-lasting illumination (73500)
- > Fiber optics provide a cool distal ring of light with no reflections or obstructions
- > Sigmoidoscopes/anoscopes can be used with a choice of illumination systems—a 6 V system with optional on/off switch and rheostat control, or with Welch Allyn's 3.5 V hand-held power sources (with adapters 73500 and 73550).
- > 73500 adapter fits most Welch Allyn 3.5 V power handles

## Sigmoidoscope/Anoscope Lighting Systems and Accessories

- |              |   |                |                                |
|--------------|---|----------------|--------------------------------|
| <b>36103</b> | Illumination System for Disposable Sigmoidoscopes (Includes: 73305, 73210, 36019 and 30200) | <b>73305</b>   | Transformer with 1.5 m cord    |
| <b>36019</b> | Fiber Optic Light Head for Disposable Sigmoidoscopes and Anoscopes                          | <b>73210</b>   | Light Handle with 0.9 m cord   |
| <b>38700</b> | Fiber Optic Light Carrier for use with Standard Anoscopes and Light Adapter                 | <b>03100-U</b> | 3.5 V Halogen Lamp (for 73550) |
| <b>73500</b> | 3.5 V Halogen Handle Adapter  | <b>30200</b>   | Insufflation Bulb, Complete    |
|              |   | <b>30130</b>   | 30 cm Suction Tube (Standard)  |
|              |   | <b>30140</b>   | 40 cm Suction Tube             |
|              |   | <b>07800-U</b> | 6 V Halogen Lamp               |

RIGID ENDOSCOPY

# Quality and Value in Rigid Endoscopy

Welch Allyn offers a wide range of rigid endoscopy products for anoscopy and sigmoidoscopy, complemented by versatile fiber optic illumination systems for optimum visualization.

The sigmoidoscope range includes both reusable and disposable instruments and a variety of sizes to accommodate differing requirements. Welch Allyn fiber optic anoscopes represent the finest in design and craftsmanship for optimal performance.



## Welch Allyn Disposable Sigmoidoscopes and Anoscopes

One of our most popular designs, the Welch Allyn Disposable Sigmoidoscope may be used in conjunction with either a 6 V illumination system or traditional Welch Allyn 3.5 V power handles.

- > Highly visible graduations to assist in positioning
- > Glare-free frosted coating ensures maximum light output at distal tip
- > Individually wrapped for safety and convenience
- > Fiber-optic illumination system provides a cool 360 degree ring of light
- > Hinged window with neoprene seal to avoid fogging during insufflation

- 53130** KleenSpec® Disposable Sigmoidoscope with Obturator (case of 100)
- 36103** Complete Illumination System for Disposables
- 73210** Light Handle with 0.9 m Cord
- 73220** Light Handle (73210) plus Transformer (73305)
- 36019** Fiber Optic Light Head
- 73305** Transformer with 1.5 m Cord
- 73500** 3.5 V Halogen Handle Adapter
- 71000-A** 3.5 V Direct Plug-in Rechargeable Handle
- 30200** Insufflation Bulb, complete
- 30210** Hydrophobic Filters (pack of 10) with Luer Lock

## Disposable Anoscope

- 53110** KleenSpec® Disposable Anoscope with Obturator (box of 25)

### Disposable Anoscope Accessories

- 73210** Light Handle with 0.9 m Cord
- 73220** Light Handle (73210) plus Transformer (73305)
- 36019** Fiber Optic Light Head
- 73305** Transformer with 1.5 m Cord



35303  
Fiber Optic Sigmoidoscope Set



## Standard and Pediatric/Stricture Sigmoidoscopes and Anoscopes

All Welch Allyn reusable sigmoidoscopes/anoscopes feature:

- > Hinged window with a neoprene seal to avoid fogging during insufflation
- > Graduated tube with removable tip for fiber cleaning (metal only)

### Standard Sigmoidoscopes

- 32810** 19 mm x 15 cm Standard Sigmoidoscope, Speculum only with Obturator
- 32820** 19 mm x 25 cm Standard Sigmoidoscope, Speculum only with Obturator
- 33830** 19 mm x 30 cm Standard Sigmoidoscope, Speculum only with Obturator
- 32830** 19 mm x 35 cm Standard Sigmoidoscope, Speculum only with Obturator
- 33220** 23 mm x 25 cm Standard Sigmoidoscope, Speculum only with Obturator

### Pediatric/Stricture Sigmoidoscopes

- 32010** 11 mm x 10 cm Fiber Optic/Stricture Speculum only with Obturator
- 32020** 11 mm x 25 cm Fiber Optic/Stricture Speculum only with Obturator
- 32410** 15 mm x 15 cm Fiber Optic/Pediatric Speculum only with Obturator
- 32420** 15 mm x 25 cm Pediatric/Stricture Speculum only with Obturator

### Rotating Anoscope

Anoscope Only 68 mm Speculum Length

- 38900** Rotating Speculum Anoscope with Obturator, 20 mm inside diameter, 22 mm outside diameter

### Fiber Optic Anoscope

Speculum Only with Obturator

- 37019** 19 mm x 7 cm Fiber Optic Anoscope, Speculum Only with Obturator
- 37023** 23 mm x 7 cm Fiber Optic Anoscope, Speculum Only with Obturator
- 37027** 27 mm x 7 cm Fiber Optic Anoscope, Speculum Only with Obturator