

# Connex<sup>®</sup> Spot Monitor Mounting and Mobility Options

Vital signs data within your reach—across care settings and workflows

Mobile stands support clinical workflows that require a shared device across multiple rooms.



	Standard Basket Stand Simple portability with a small footprint	Work Surface and Accessory Management Stand Organized storage and cable management, plus a compact mobile work surface footprint	Accessory Power Management Stand Mobile work surface and integrated battery to power through a busy day
--	--	---	--

<b>Part Number*</b>	7000-MS3	7000-MWS	7000-APM
<b>Height</b>	49" (124.46 cm) total 38" (96.52 cm) to top of basket handle	47" (119.38 cm) total 39.25" (99.7 cm) to work surface	47" (119.38 cm) total 41" (104.14 cm) to work surface
<b>Base Diameter</b>	23" (58.42 cm)	23" (58.42 cm)	23" (58.42 cm)
<b>Work Surface Area</b>	N/A	11.62" (29.51 cm) wide 11.25" (28.56 cm) deep	11.5" (29.21 cm) wide 10" (25.4 cm) deep
<b>Accessory Basket</b>	12.5" (31.75 cm) wide 8.25" (20.96 cm) deep 7.25" (18.42 cm) high	11.75" (29.85 cm) wide 8.25" (20.96 cm) deep 7.5" high (19.05 cm) 4 accessory bays	11.75" (29.85 cm) wide 8.25" (20.96 cm) deep 7.5" high (19.05 cm) 4 accessory bays
<b>Weight</b>	13.8 kg / 30.5 lb (fully loaded with device and accessories)	15.4 kg / 34 lb (fully loaded with device and accessories)	18.4 kg / 40 lb (fully loaded with device and accessories)
<b>Power Accessory Functionality</b>	N/A	N/A	9 cell Lithium Ion Battery (Field replaceable, integrated power supply, integrated work light, battery life indicator display, two USB ports)

\*Connex Spot Monitor sold separately.

Wall mount options are ideal for care settings that need a dedicated device in every exam room, while desktop stands can be moved as needed.



#### Desktop Stand

Compact and portable tabletop mount suitable for placement on any work surface

#### Wall Mount with Basket

Easy-to-access storage and cable management

#### Integrated Wall System

Complete set of physical exam tools in one integrated package

	Desktop Stand	Wall Mount with Basket	Integrated Wall System
<b>Part Number(s)*</b>	7000-DST	7000-GCX	77794-MCSM**, 77794-2MCSM***
<b>Dimensions</b>	Height 9.5" (24.13 cm) Width 7.25" (18.42 cm) Depth 4.75" (12.07 cm)	Wall Channel: 19" (48.26 cm) Pivot Arm: 8" (20.32 cm) Wall Mount Utility Basket: 12.75" (32.39 cm) wide x 7.5" (19.05 cm) deep x 5.9" (14.99 cm) high	Height: 12" (30.48 cm) Length: 34" (86.36 cm) Depth*: 2" (5.08 cm) *Depth dimension is for wall board only

The Connex Spot Monitor uses a standard VESA mount. Aftermarket mounting options are available.



Learn more at [www.welchallyn.com](http://www.welchallyn.com).

\* Connex Spot Monitor sold separately.

\*\* 77794-MCSM model includes mounting bracket for Connex Spot Monitor, Coaxial Ophthalmoscope (#11720), MacroView™ Otoscope (#23810), KleenSpec® Disposable Specula dispenser, vital signs accessory basket, vital signs monitor power supply (#7000-PS) and wall mount panel.

\*\*\* 77794-2MCSM model includes mounting bracket for Connex Spot Monitor, PanOptic™ Ophthalmoscope with cobalt-blue filter and corneal viewing lens (#11820), MacroView Otoscope (#23810), KleenSpec Disposable Specula dispenser, vital signs accessory basket, vital signs monitor power supply (#7000-PS) and wall mount panel.